OPTIC & SCIENTIFIC TOYS

www.navir.it

ITEM 8050

03-MS082

ZOOMSCOPE

microscope

A journey into the microcosmos

Made in Italy

NAVIR snc

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WARNING ! : This functions with battery type LR1 (size N) , 1,5 V.

Only a battery of the same type or equivalent to that recommended should be used. The rundown battery must be removed and thrown in the special differentiated waste container.

Read and keep for future reference

- Never recharge non-rechargeable batteries.
- Rechargeable batteries must only be recharged under adult supervision.
- Rechargeable batteries must be removed from the toy before being recharged.
- The charging clips must not be short-circuited.

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max 1,2 V - 0,22 A

Brief historical outline




Lenses were already known to the Greeks and Romans who used them to light fires (burning lenses). It was in 1300 that lenses were first produced in Italy to manufacture eyeglasses, a production that developed greatly in Holland.

It was precisely in Holland that in 1600 it was discovered that, by adequately combining together different lenses, it was possible to magnify either very distant or very small objects. Thus were the microscope and the telescope born, two instruments that have enabled humankind to progress in the knowledge of the world we live in.

The microscope, first manufactured with a single lens

Left: Burning lens. Bottom: Galileo Galilei and a Dutch eyeglasses merchant (16th century).

Left: Antony van Leeuwenhoek simple microscope. Right: Robert Hooke's compound microscope. Bottom: Hooke's drawing representing the cells in the bark of cork oak.

(simple microscope) and later with two (compound microscope), has enabled us to study the hidden world which is not visible to the naked eye and, in particular, made finally visible formerly unknown forms of life.




Among the first to use the microscope were the Italian scientist Galileo Galilei (1564-1642) who actually gave the instrument its name (from the Greek mikros = little and skopein = to see), the Dutch Antony van Leeuwenhoek (1632-1723) and the English Robert Hooke (1635-1703). The latter discovered the existence of cells by observing a piece of cork oak bark.

Although Zoomscope may seem just a toy, it should be remembered that many important discoveries were made with instruments with like power and characteristics.

The Zoomscope



The Zoomscope consists of an eye lens, an eyepiece and a mobile lens that allows to vary the degree of magnifying power. The magnifying power can be varied by turning a wheel situated on the side of the instrument.

Usually one starts with the lowest magnification to increase it gradually, until an optimal view of the object is achieved.


ENGLISH - 1


Top: The lighting system is situated in its housing.
Above: light switching on and off.
Right: battery change.

Like all microscopes, Zoomscope produces an upside-down image of the observed object.

Specimen lighting


The Zoomscope allows two kinds of observation. If the material is transparent, it can be observed in transmitted light, i.e. with the light that passes through the material. This kind of viewing is achieved by aiming the

instrument towards a light source such as the sky or a lamp (1).

Vice versa, if the material is opaque, one should use incident light. In this case the built-in lighting system of the instrument should be switched on, thus lighting whatever has been placed within the field of vision (2). When the instrument lighting system is not in use, it can be stored in its proper housing.

The battery should be inserted as shown in the figure. To switch the bright spotlight from on to off, turn the small handle on the side opposite to the light.


Focusing

Focusing is achieved by slowly rotating the transparent part (or, alternatively, by rotating the instrument). If the object you are viewing is flat, it can be brought into

Caution!

Never aim the Zoomscope (or any other optical instrument) towards the sun, for it can cause eye damage.

focus with just one operation. Vice versa, if the object is thick, only one section at a time can be focused. Therefore, to see other sections of the object that are not on the same plane, the focus must be adjusted continuously.


Slides

Sometimes the specimens are better observed on a slide or between two slides. Zoomscope is supplied with two slides.

To observe a slide, the stage must be assembled to the instrument. The specimen is then applied to the slide (either dry or moistened with a drop of water). Finally, the slide is covered with another slide and the two are fixed together with adhesive tape.

The material can be better observed if it has previously been

Below, pictures of a Clostridium, focused in three different ways.

stained, for instance with iodine tincture.

The micrometer

One of the two supplied slides presents two small scales to measure the size of the observed object. One has a notch every

half millimeter and the other has one every 32nd of an inch.

To measure the observed object, just place it next to the scale, which is used as if it were a ruler.


Viewing


The Zoomscope allows the observation of any object, live or inanimate, placed within the field of vision and brought into focus. Hereunder we are listing a number of things that can be interesting to observe:

Insects

Insects offer many opportunities for observation. The one at the left is a small fly.

Onion skin.

Leaves and other vegetable parts

Vegetables should be sliced as thin as possible. To this purpose professionals use an instrument called a microtome. You can use a sharp knife or a razor blade (always with the assistance of an adult).

The picture shows the thin skin of an onion.

Crystals

Dissolve a small amount of kitchen salt in little boiling water and stir for a long time. In this way you obtain a saturated salt solution. Apply one or two drops of this solution to the slide and allow all the water to evaporate. What is left after all the water has evaporated are the typical cubic crystals of sodium chloride, i.e. of common salt. You can follow the same procedure with sugar or any other crystalline substance.

Hair

Observe the difference between straight, curly or frizzy hair.

Stamps, notes, documents

The Zoomscope can be used to study and verify the authenticity of every kind of printed matter and documents.

Fingerprints

Police all over the world use the microscope to study and compare fingerprints.

Fabrics

Fabric weaves can be observed in detail.

Printed matter

Color prints are an interesting object to observe with the Zoomscope.

With it you can learn how an infinite range of colors can be obtained simply by different combinations of the four basic colors: cyan, magenta, yellow and black.

Above, Daphnia longispina.

Sodium chloride (kitchen salt) crystals.

From the left: bank note and stamp details.

cyan

magenta

yellow

black

4 colors

Useful accessories

In many toyshops or shops of naturalists' items you can find many useful articles, such as:

- slides, cover slips and other instruments for the preparation of specimens;
- ready-made slides with specimens of small animals, animal and vegetable tissue etc.;
- manuals illustrating the various techniques used in the preparation of specimens.

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